




**Laboratory**  
Design & Supply

## Counter Top Double Bowl Sink



This unit is designed with marine edge on horizontal surfaces to avoid countertop spilling and a backsplash unit. Uniquely designed one piece molding with integral sink outlets, edge and backsplash. Unit has attached 1-1/8" plywood for easy affixing to any counter or cabinet top unit. Standard color is non-glare black.

SYMBOL	A	B	C	D	E	F	G	J	K	L	M	N
W31063-000	18"	15"	10"	49"	30"	3"	2.6875"	8.1875"	5"	8.1875"	4"	1.375"