




Laboratory Design & Supply

SINGLE BOWL CENTER DRAIN POLYOLEFIN SINKS

Single Bowl Center Drain

This series of sinks allows for uniform drainage from all locations on the sink. Colors available are Black, White, and Natural.



Part Number	A	B	C	D	E	F	G	H	OUTLET
W31415-000	11	11	7	13-1/2	13-1/2	5-1/2	9-3/4	9-3/4	1-1/2
W31418-000	12	8	8	14-1/2	10-1/2	4	6	11	1-1/2
W31440-000	16	6	6	18	8	2-3/4	7-3/4	8-1/2	1-1/2
W31449-000	16	6	6	18-1/2	8-1/2	2-3/4	7-1/2	8-1/2	1-1/2
W31466-000	18	15	12	20-1/2	17-1/2	7-1/2	9	15	1-1/2
W31473-000	21	18	12	23-1/2	20-1/2	9	10-1/2	15	1-1/2
W31477-000	23	17	8	25-1/2	19-1/2	9	10-1/2	11	1-1/2
W31478-000	23	17	11-1/2	25-1/2	19-1/2	8-3/4	11-1/2	11	1-1/2
W31481-000	24	20	10	26-1/2	22-1/2	10	12	13	1-1/2
W31483-000	24	20	12	26-1/2	22-1/2	10	12	15	1-1/2